## Choose your GO<sub>2</sub> Accessories



#### GO2CC

Black neoprene carrying case with belt loop



GO2L

Black 20" Lanyard



GO2R

Retractable Clip Holder



## A Message to Oxygen Users

## With Knowledge Comes Peace of Mind

Having the necessary information is essential to successfully managing your oxygen to stay within your clinician-prescribed oxygen saturation and heart rate range. With **GO<sub>2</sub>**, you're in control:

- Adjust your oxygen to meet your daily activity levels.
- Get instant feedback on the effectiveness of your breathing techniques.
- Improve communication with your clinician using trusted Nonin Technology.

"I won't leave home without it." -Edna Fiore, COPD Patient/Advocate

For oxygen patients who love to travel, the new FAA rules allow the use of Portable Oxygen Concentrators (POCs) on commercial airlines.

## Go Further! Independence Restored.

A chronic lung disease diagnosis often requires that adjustments be made in your daily activities. However, it does not have to mean a life of **limitation!** Whether you're looking for peace of mind, to go out for a night at the movies, or to up your activity level by walking those extra blocks around the neighborhood — the GO<sub>2</sub> fingertip oximeter is the perfect companion! Wherever you decide to go, GO<sub>2</sub> tucks away in your purse or pocket and provides the assurance that you are within your prescribed oxygen saturation and heart rate range — anywhere, anytime.

> An oximeter gives the patient immediate feedback, a lot of reassurance — reduces anxiety. –Patty Koff, RRT

#### **Specifications**

0% to 100 Oxygen Saturation Display Range Pulse Rate Display Range 18 to 321 pulses per minute

Oxygen Saturation Declared Accuracy

Range (Arms\*) Low Perfusion Oxygen Saturation Declared

Accuracy Range (Arms\*)

Pulse Rate Declared Accuracy Range(Arms\*)

Low Perfusion Pulse Rate Declared Accuracy

Range (Arms\*)

Storage/Transportation

Battery Life (Storage)

\*± 1 Arms represents approximately 68% of measurements.

+41° to +104°F (5°C to +40°C) Temperature (Operating) Storage/Transportation -22° to +158°F (-30° to +70°C)

Device temperature will not exceed 41°C as measured in a controlled environment.

**Humidity** (Operating) 10% to 90% relative humidity,

non-condensing

non-condensing

10% to 95% relative humidity,

70-100% SpO2 ± 2 digits

70-100% SpO2 ± 2 digits

20-250 BPM ± 3 digits

40-240 ± 3 digits

Operating Altitude Up to 40,000 feet / 12.192 meters Battery Life (Continuous): Approximately 2,400 spot

checks based on ~21 hours of

operation using one AAA-size alkaline battery, calculated at 30

seconds per use.

6 months minimum

2 year warranty. Backlit LCD for low light viewing. Intended for Home Health Care Only.

Specifications are subject to change without notice.



#### Nonin Medical, Inc.

13700 1st Avenue North • Plymouth MN • 55441-5443 • USA Phone: +1.763.553.9968 • Toll Free: 877.577.2635 • Fax: +1.763.577.5521 E-mail: info@nonin.com



## www.go2nonin.com









All trademarks are the property of Nonin Medical. Inc unless otherwise noted







## **Empower Yourself**

Nonin Medical, the inventor of fingertip pulse oximetry, introduces the GO<sub>2</sub> fingertip pulse oximeter. The affordable GO<sub>2</sub> empowers you with the knowledge of reliable oxygen saturation and heart rate values without compromising the quality and accuracy that have made Nonin the most trusted name in fingertip pulse oximetry.

Live life on the GO!



### Go Ahead and Get Active!

A comprehensive fitness program is essential to maintaining lung health and quality of life if you have been diagnosed with a chronic lung disease. The highly durable  $GO_2$  fingertip oximeter is the perfect oxygen saturation and heart rate monitoring device to complement your clinician-devised exercise program. *Know when to increase your liter flow in an instant!* Compact and easy-to-use, simply place  $GO_2$  on your finger for accurate oxygen saturation and heart rate feedback. Go further with  $GO_2$ !



Patient/Advocate

## **Quality and Support DO Matter**

With so many throwaway products available these days, it's important to choose a product that not only is consistently accurate, but is also built to last. At Nonin, we are dedicated to providing the highest quality products available at a great value.

## **Superior Service and Support**

Nonin will be here long after the sale and our knowledgeable customer service representatives are always available to answer your questions. Patients should all have oximeters to guide their oxygen use during all activities of daily living.

-Dr. Thomas L. Petty, MD

## A Message to Clinicians

## **Trusted Accuracy**

Accurate reading of  $SpO_2$  values is among the most important factors for successful management of LTOT care. With the  $GO_2$  fingertip oximeter, engineered with Nonin Medical's proven PureSAT®  $SpO_2$  technology (the same technology that's in our medical grade oximeters), you can always trust the oxygen saturation and heart rate values reported by your patients. With over 22 years of manufacturing, engineering and research experience, we have a long-standing reputation for excellence. In fact, 96% of pulmonary rehab facilities use Nonin.

## **Partnering with your Patients**

Having a fingertip oximeter that is always accurate helps patients stay within their prescribed oxygen saturation and heart rate levels — enabling better management of their care by giving them the information they need to adjust their oxygen to meet their activity levels. With the power and confidence to manage their oxygen, the  $GO_2$  fingertip oximeter empowers patients to respond in situations where they find themselves not feeling 100%. Visit www.go2nonin.com to learn more about how the  $GO_2$  can benefit your patients.

Every COPD patient should have their own oximeter. I treat it just the same as I would having their own blood pressure cuff.

—Brenda Crowe, RT

## **Promoting Exercise Compliance**

With the reassurance of the  $GO_2$  fingertip oximeter by their side, patients will be more apt to comply with an exercise regimen.

# **Empowering Patients to Take Control** of their Lives

A dependable and accurate tool like the **GO**<sub>2</sub> fingertip oximeter will empower patients with the knowledge and confidence that they can manage this disease and still live a fulfilling life.

## A Message for All

Not All Oximeters are Created Equal!

- Nonin welcomes scrutiny of our independent testing results and they are available for all to see!
  Check out why Nonin is the most trusted name in fingertip oximetry at www.go2nonin.com.
- Anyone can make low perfusion claims, but ONLY Nonin has claims for fingertip pulse oximeters cleared by the FDA. Check out our accuracy specifications at www.go2nonin.com. Demand the same from the imitators.
- Nonin's quality standards require every **GO**<sub>2</sub> be able to withstand 25 drops from one meter onto a hard surface. Try it with a **GO**<sub>2</sub> imitator and be prepared to pick up the pieces!
- Home environments can get messy **GO**<sub>2</sub> exceeds IP33 water ingress testing.
- The GO<sub>2</sub> is NOT a throw-away product. If something happens, it can be repaired and like new in record time. Nonin is a company that will be here tomorrow.
- At Nonin, you'll speak with a real person right away. The voice you hear represents the pride in our products. We go above and beyond to have satisfied customers.